



This is to certify that units of the lot number below were tested and found to comply with the specifications of the grade listed. Products are tested to meet the analytical requirements of the noted grade. The following information is the actual analytical results obtained, for reference.

Part Number:	CL-F075P		
Lot.No.:	202312010101		
Specification:	Size: <input type="checkbox"/> 25cm ² <input checked="" type="checkbox"/> 75cm ² <input type="checkbox"/> 175cm ² <input type="checkbox"/> 225cm ²		
	Surface treatment: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Cap type: <input checked="" type="checkbox"/> Plug cap <input type="checkbox"/> Vent cap		
	Bottle Material: <input checked="" type="checkbox"/> PS		Cap Material: <input checked="" type="checkbox"/> HDPE
Product Name:	Cell Culture Flask		
Package:	5pcs/bag, 20bags/case		
Shelf Life:	5 years		
Manufacturer:	PakGent Bioscience (Suzhou) Co., Ltd.		
Inspection Date:	21st, May, 2026	Report Date	21st, May, 2026
Inspection items and result:			
Cosmetic:	No black spots, no cuts, clear and transparent		
Package:	Accurate quantity and label, inner package without damage		
Dimension:	Planeness Bottom ≤ 0.5mm, Cover ≤ 0.5mm		
Function test:	N/A		
Endotoxin:	<0.005EU/ml		
RNase/DNase:	DNase ≤ 4x10 ⁻⁶ Kunitz units/μl; RNase ≤ 10 ⁻⁹ Kunitz units/μl		
ATP:	2x10 ⁻¹² mg/μl		
PCR Inhibitors:	Not detectable		
Human DNA:	<30pg		
Sterile:	Irradiation, SAL 10 ⁻⁶		
Conclusion: Comply with Specification			

This signature indicates that the above test has met all quality specification and has been reviewed by a quality representative.

Title: QA Manager

Signature: Alvin Cao

PakGent Bioscience (Suzhou) Co., Ltd.

Tel: +86 (0) 512 6670 0017

Address: Floor 1-3, Building 33, Alloy Industrial Park, Wuzhong District, Suzhou, China.

PakGent Bioscience (Shanghai) Co., Ltd.

Tel: 400-167-1717

Address: Building C, No 888, Lake west road, Lin gang section, Free Trade Zone, Shanghai, China.

PakGent Bioscience (Changshu) Co., Ltd.

Tel: +86 (0) 512 5298 8330

Address: No.9 Lecheng Road, Bixi Street, Changshu, Jiangsu Province, China.