


Technical Datasheet

Article - No:	TDS-PCRP-20IS	Revision: A0
Remarks	single use only	
Specification		
Product name	PCRP-20IS,0.2ml PCR 96-well plate	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Nominal Capacity	0.2ml	
With skirt	sub semi-skirt or high-skirt	
Material Characteristics	polypropylene (PP), premium surface	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per box: 10	
	box per case: 10	
Storage	dry condition:40%-60% rel. humidity	
Approved by PakGent quality department		