


Technical Datasheet

Article - No:	TDS-EFC-250P	Revision: A0
Remarks	single use only	
Specification		
Product name	EFC-250P,250ml,Erlenmeyer Flask,Plug Cap,Flat bottom	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	250ml	
Material Characteristics	bottle:Polycarbonate(PC)	
	cap:High Density Polyethylene(HDPE)	
Sterile	Yes	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	60 months at controlled storage conditions	
Packaging	pcs per case: 12	
Storage	Storage Temperature:20°C-25°C,humidity40%-60%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		