


## Technical Datasheet

|   |  |              |
|---|--|--------------|
| <b>Article - No:</b>                          | TDS-CL-FT175V  | Revision: A0 |
| <b>Remarks</b>                                | single use only  |              |
| Specification                                 |  |              |
| <b>Product name</b>                           | CL-FT175V, 175cm <sup>2</sup> , Cell Culture Flask, TC treated, Vent cap           |              |
| <b>Production conditions</b>                  | Class 10,000 clean room production   |              |
|   | No release agent is used   |              |
| <b>Classification</b>                         | CE/IVD   |              |
| <b>Dimensions</b>                             |  |              |
| <b>Colour</b>                                 | Clear  |              |
| <b>Graduation</b>                             | moulded rings  |              |
| <b>Nominal Capacity</b>                       | 175cm <sup>2</sup>   |              |
| <b>Material Characteristics</b>               | bottle: polystyrene (PS)   |              |
|   | cap: High Density Polyethylene (HDPE)  |              |
| <b>Sterile</b>                                | Yes  |              |
| <b>Product Characteristics</b>                | biologically inert   |              |
|   | no heavy metals according to EC/94/62  |              |
| <b>Purity Grade</b>                           | free of RNase and DNase  |              |
|   | free of pyrogen  |              |
|   | free of PCR inhibitors   |              |
|   | free of human DNA  |              |
|   | free of endotoxins   |              |
|   | free of ATP  |              |
| <b>Shelf - Life</b>                           | 60 months at controlled storage conditions   |              |
| <b>Packaging</b>                              | pcs per bag: 5   |              |
|   | bag per case: 6  |              |
| <b>Storage</b>                                | Storage Temperature: 20°C-25°C, humidity 40%-60%                                   |              |
|   | Keep clean, ventilated and out of direct sunlight                                  |              |
| <b>Approved by PakGent quality department</b> |  |              |