


Technical Datasheet

Article - No:	TDS-WMPB1000A	Revision: A0																									
Remarks	single use only																										
Specification																											
Product name	WMPB1000A, 1000ml Packaging Bottle, HDPE																										
Production conditions	Class 100,000 clean room production																										
	No release agent is used																										
Classification	CE/IVDD																										
Dimensions																											
	<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th>名称</th> <th>内容</th> <th>瓶身材料</th> <th>瓶身外径 D1 (mm)</th> <th>瓶身高度 H1 (mm)</th> <th>瓶口内径 D2 (mm)</th> <th>瓶盖外径 D3 (mm)</th> <th>瓶盖高度 H2 (mm)</th> <th>产品总高 连盖 (mm)</th> </tr> </thead> <tbody> <tr> <td rowspan="2">1000ml 试剂瓶</td> <td>PP</td> <td>PP</td> <td>φ91.9±0.5</td> <td>193.9±0.5</td> <td>φ43.7±0.4</td> <td>φ57.5±0.4</td> <td>16.8±0.4</td> <td>196.0±1.0</td> </tr> <tr> <td>PE</td> <td>PE</td> <td>φ92.2±0.5</td> <td>193.9±0.5</td> <td>φ43.7±0.4</td> <td>φ57.5±0.4</td> <td>16.8±0.4</td> <td>196.0±1.0</td> </tr> </tbody> </table>		名称	内容	瓶身材料	瓶身外径 D1 (mm)	瓶身高度 H1 (mm)	瓶口内径 D2 (mm)	瓶盖外径 D3 (mm)	瓶盖高度 H2 (mm)	产品总高 连盖 (mm)	1000ml 试剂瓶	PP	PP	φ91.9±0.5	193.9±0.5	φ43.7±0.4	φ57.5±0.4	16.8±0.4	196.0±1.0	PE	PE	φ92.2±0.5	193.9±0.5	φ43.7±0.4	φ57.5±0.4	16.8±0.4
名称	内容	瓶身材料	瓶身外径 D1 (mm)	瓶身高度 H1 (mm)	瓶口内径 D2 (mm)	瓶盖外径 D3 (mm)	瓶盖高度 H2 (mm)	产品总高 连盖 (mm)																			
1000ml 试剂瓶	PP	PP	φ91.9±0.5	193.9±0.5	φ43.7±0.4	φ57.5±0.4	16.8±0.4	196.0±1.0																			
	PE	PE	φ92.2±0.5	193.9±0.5	φ43.7±0.4	φ57.5±0.4	16.8±0.4	196.0±1.0																			
Colour	Amber																										
Nominal Capacity	1000ml																										
Material Characteristics	bottle: High Density Polyethylene (HDPE)																										
	cap: polypropylene (PP)																										
Performance parameters	Negative pressure: -70kpa ~ -90kpa																										
	Ethanol loss rate: ≤ 0.5%																										
	water loss rate: ≤ 0.2%																										
Purity Grade	free of RNase and DNase																										
	free of pyrogen																										
	free of PCR inhibitors																										
	free of human DNA																										
	free of endotoxins																										
Shelf - Life	60 months at controlled storage conditions																										
Packaging	pcs per case: 50																										
Storage	dry condition: 40%-60% rel. humidity																										
Approved by PakGent quality department																											