

Technical Datasheet

Article - No:	TDS-STSS-509CR	Revision: A0
Remarks	single use only	
Specification		
Product name	STSS-509CR,0.5ml Screw Cap Tubes	
Production conditions	Class 10,000 clean room production	
	No release agent is used	
Classification	CE/IVDD	
Dimensions	<p>The drawing shows a side view of the tube with a screw cap. Dimensions are: total length 45.00, length to the start of the cap 37.50, diameter of the tube body 12.50, and diameter of the cap 9.80. A top view of the cap shows an outer diameter of 12.50.</p>	
Colour	Clear tube,Red cap	
Nominal Capacity	0.5ml	
Material Characteristics	polypropylene (PP)	
Performance parameters	Storage Temperature: -80°C~121°C	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	60 months at controlled storage conditions	
Packaging	pcs per bag:500	
	bag per case:10	
Storage	dry condition:40%-60% rel. humidity	
Approved by PakGent quality department		