


Technical Datasheet

Article - No:	TDS-STSS-511CNS	Revision: A0
Remarks	single use only	
Specification		
Product name	STSS-511CNS,2.0ml Screw Cap Tubes,individual pack	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVDD	
Dimensions		
Colour	Clear	
Nominal Capacity	2.0ml	
Material Characteristics	polypropylene (PP)	
Sterile	Yes	
Performance parameters	Storage Temperature: -80°C~121°C	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	60 months at controlled storage conditions	
Packaging	pcs per bag:500	
	bag per case:10	
Storage	dry condition:40%-60% rel. humidity	
Approved by PakGent quality department		