

# Large volume tips

## Guidelines For Use

No. : PGLVT20230418E



PakGent tips are made of 100% pure virgin medical grade polypropylene, so you will not find recycled plastics in our tips. We pre-test the resin to ensure no extractables are present.

Every lot of tips is lab tested and certified to be free of detectable Nuclease (Rnase/Dnase), Adenosine Tri-phosphate, and Endotoxin (pyrogen) contamination.

Large volume (5 ml, 10 ml) pipette tips can dispense large volume samples with the same accuracy as smaller volume tips.

These can reduce or eliminate the use of fragile glass and plastic serological pipettes. Hands-free tip loading is now possible when using 5 ml and 10 ml pipettes.

Product characteristics	Description
Material	Polypropylene.
Sterilization grade	SAL 10 <sup>6</sup> .
Certified quality	Our products are made to ISO standards and we are audited lot to lot to ensure you continuous quality.
Lab tested	PakGent tips are developed using our own in-house laboratory and tested on real pipettors. We use actual lab results on real pipettes.
Packing	Available of bag and rack pack.
Mainstream brand pipettes	Fit Eppendorf, Biohit, DLAB, Thermo, Brand Pipette etc.

### 5ml, clear, fit Eppendorf, Biohit, DLAB Pipette




Part NO.	With Filter	Sterile	Pack
PT-5000R-E	✗	✓	24/rack, 30racks/case
PT-5000RF-E	✓	✓	24/rack, 30racks/case
PT-5000B-E	✗	✗	100/bag, 20bags/case
PT-5000BF-E	✓	✗	100/bag, 20bags/case

### 10ml, clear, fit Eppendorf, Biohit, DLAB Pipette



Part NO.	With Filter	Sterile	Pack
PT-10000R-E	✗	✓	24/rack, 40racks/case
PT-10000RF-E	✓	✓	24/rack, 40racks/case
PT-10000B-E	✗	✗	100/bag, 10bags/case
PT-10000BF-E	✓	✗	100/bag, 10bags/case

### 5ml, clear, fit Thermo, Brand Pipette



Part NO.	With Filter	Sterile	Pack
PT-5000R-T	✗	✓	40/rack, 40racks/case
PT-5000RF-T	✓	✓	40/rack, 40racks/case
PT-5000B-T	✗	✗	100/bag, 20bags/case
PT-5000BF-T	✓	✗	100/bag, 20bags/case

### 10ml, clear, fit Thermo, Brand Pipette



Part NO.	With Filter	Sterile	Pack
PT-10000R-T	✗	✓	24/rack, 40racks/case
PT-10000RF-T	✓	✓	24/rack, 40racks/case
PT-10000B-T	✗	✗	100/bag, 10bags/case
PT-10000BF-T	✓	✗	100/bag, 10bags/case

## Compatibility List

PIPETTE BRAND		Eppendorf		Thermo		Biohit		DLAB		Brand	
		Research plus	Research plus	Finnpipette F3	Finnpipette F3	Proline Plus	Proline Plus	HiPette	HiPette	Transferpette® S	Transferpette® S
PRODUCT NO.	PIPETTE VOLUME	0.5-5ml	1-10ml	0.5-5ml	1-10ml	500-5000ul	1-10ml	1000-5000ul	2-10ml	500-5000ul	1000-10000ul
PT-5000B-E	500-5000ul	√				√		√			
PT-5000BF-E	500-5000ul	√				√		√			
PT-5000R-E	500-5000ul	√				√		√			
PT-5000RF-E	500-5000ul	√				√		√			
PT-10000B-E	1000-10000ul		√				√		√		
PT-10000BF-E	1000-10000ul		√				√		√		
PT-10000R-E	1000-10000ul		√				√		√		
PT-10000RF-E	1000-10000ul		√				√		√		
PT-5000B-T	500-5000ul			√						√	
PT-5000BF-T	500-5000ul			√						√	
PT-5000R-T	500-5000ul			√						√	
PT-5000RF-T	500-5000ul			√						√	
PT-10000B-T	1000-10000ul				√						√
PT-10000BF-T	1000-10000ul				√						√
PT-10000R-T	1000-10000ul				√						√
PT-10000RF-T	1000-10000ul				√						√

**PakGent®**

PakGent Bioscience (Suzhou) Co., Ltd.

Address: No. 9 Lecheng Road, Bixi Street, Changshu, Jiangsu Province, China

Tel: +86 (0) 512 6670 0017

Email: [info@pakgentbio.com](mailto:info@pakgentbio.com)

Website: [www.pakgentbio.com](http://www.pakgentbio.com)



Catalog