


Technical Datasheet

Article - No:	TDS-NMPB125P	Revision: A0																							
Remarks	single use only																								
Specification																									
Product name	NMPB125P, 125ml Packaging Bottle, PP																								
Production conditions	Class 100,000 clean room production																								
	No release agent is used																								
Classification	CE/IVDD																								
Dimensions	 <table border="1" style="margin-left: auto; margin-right: auto;"> <thead> <tr> <th>location name</th> <th>material</th> <th>O.D. D1 (mm)</th> <th>Height H1 (mm)</th> <th>I.D. D2 (mm)</th> <th>Cap O.D. D3 (mm)</th> <th>Cap H H2 (mm)</th> <th>Total H (mm)</th> </tr> </thead> <tbody> <tr> <td rowspan="2">125ml</td> <td>PP</td> <td>φ50.4</td> <td>97.9</td> <td>φ17.5</td> <td>φ28.4</td> <td>15.1</td> <td>99.1</td> </tr> <tr> <td>PE</td> <td>φ50.4</td> <td>97.9</td> <td>φ17.5</td> <td>φ28.4</td> <td>15.1</td> <td>99.1</td> </tr> </tbody> </table>		location name	material	O.D. D1 (mm)	Height H1 (mm)	I.D. D2 (mm)	Cap O.D. D3 (mm)	Cap H H2 (mm)	Total H (mm)	125ml	PP	φ50.4	97.9	φ17.5	φ28.4	15.1	99.1	PE	φ50.4	97.9	φ17.5	φ28.4	15.1	99.1
	location name	material	O.D. D1 (mm)	Height H1 (mm)	I.D. D2 (mm)	Cap O.D. D3 (mm)	Cap H H2 (mm)	Total H (mm)																	
125ml	PP	φ50.4	97.9	φ17.5	φ28.4	15.1	99.1																		
	PE	φ50.4	97.9	φ17.5	φ28.4	15.1	99.1																		
Colour	Clear																								
Nominal Capacity	125ml																								
Material Characteristics	bottle: polypropylene (PP)																								
	cap: High Density Polyethylene (HDPE)																								
Performance parameters	Negative pressure: -70kpa ~ -95kpa																								
	Ethanol loss rate: ≤ 0.5%																								
	water loss rate: ≤ 0.2%																								
Purity Grade	free of RNase and DNase																								
	free of pyrogen																								
	free of PCR inhibitors																								
	free of human DNA																								
	free of endotoxins																								
Shelf - Life	60 months at controlled storage conditions																								
Packaging	pcs per case: 500																								
Storage	dry condition: 40%-60% rel. humidity																								
Approved by PakGent quality department																									