

Technical Datasheet

Article - No:	TDS-SPTB05	Revision: A0
Remarks	single use only	
Specification		
Product name	SPTB05,5ML Serological Pipette,Bulk	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	<p>The drawing shows a 5ml serological pipette with a total length of 335±2 mm. The tip has a diameter of 30±0.5 mm and a length of 2.0±0.5 mm. The main body has a diameter of 6.70±0.1 mm. The graduation markings are printed on the body, with major markings every 1 ml and minor markings every 0.02 ml. The pipette is labeled 'Serological' and '5ml'.</p>	
Colour	Clear	
Graduation	Printing	
Nominal Capacity	5ML	
Sterile	Yes	
Material Characteristics	Polystyrene (PS), premium surface	
	filter:high density polyethylene (PE)	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	60 months at controlled storage conditions	
Packaging	pcs per bag: 50	
	bag per case: 10	
Storage	Storage Temperature:2°C-40°C, humidity40%-80%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		