


Technical Datasheet

| | | |
|---|---|--------------|
| Article - No: | TDS-SPTS02 | Revision: A0 |
| Remarks | single use only | |
| Specification | | |
| Product name | SPTS02,2ML Serological Pipette,Individual packaging | |
| Production conditions | Class 100,000 clean room production | |
| | No release agent is used | |
| | Motor molding | |
| Classification | CE/IVD | |
| Dimensions |  | |
| Colour | Clear | |
| Graduation | Printing | |
| Nominal Capacity | 2ML | |
| Sterile | Yes | |
| Material Characteristics | Polystyrene (PS), premium surface | |
| | filter:high density polyethylene (PE) | |
| Product Characteristics | biologically inert | |
| | no heavy metals according to EC/94/62 | |
| Purity Grade | free of RNase and DNase | |
| | free of pyrogen | |
| | free of PCR inhibitors | |
| | free of human DNA | |
| | free of endotoxins | |
| | free of ATP | |
| Shelf - Life | 60 months at controlled storage conditions | |
| Packaging | pcs per bag: 50 | |
| | bag per case: 20 | |
| Storage | Storage Temperature:2°C-40°C,humidity40%-80% | |
| | Keep clean,ventilated and out of direct sunlight | |
| Approved by PakGent quality department | | |