



## Technical Datasheet

<b>Article - No:</b>	TDS-STSS-511AM	Revision: A0
<b>Remarks</b>	single use only	
<b>Specification</b>		
<b>Product name</b>	STSS-511AM,2.0ml Screw Cap Tubes	
<b>Production conditions</b>	Class 100,000 clean room production	
	No release agent is used	
<b>Classification</b>	CE/IVDD	
<b>Dimensions</b>	 	
<b>Colour</b>	Amber tube,Amber cap	
<b>Nominal Capacity</b>	2.0ml	
<b>Material Characteristics</b>	polypropylene (PP)	
<b>Performance parameters</b>	Storage Temperature: -80°C~121°C	
<b>Product Characteristics</b>	biologically inert	
	no heavy metals according to EC/94/62	
<b>Purity Grade</b>	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
<b>Shelf - Life</b>	60 months at controlled storage conditions	
	60 months at controlled storage conditions	
<b>Packaging</b>	pcs per bag:500	
	bag per case:10	
<b>Storage</b>	dry condition:40%-60% rel. humidity	
<b>Approved by PakGent quality department</b>		