


Technical Datasheet

Article - No:	TDS-WMPB060	Revision: A0																						
Remarks	single use only																							
Specification																								
Product name	WMPB060HA,60ml Packaging Bottle,HDPE, Amber																							
Production conditions	Class 100,000 clean room production																							
	No release agent is used																							
Classification	CE/IVDD																							
Dimensions																								
	<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th>location name</th> <th>material</th> <th>O. D. D1 (mm)</th> <th>Height H1 (mm)</th> <th>I. D. D2 (mm)</th> <th>Cap O.D D3 (mm)</th> <th>Cap H H2 (mm)</th> <th>total H (mm)</th> </tr> </thead> <tbody> <tr> <td rowspan="2">60ml</td> <td>PP</td> <td>$\phi 39.0 \pm 0.3$</td> <td>81.5 ± 0.3</td> <td>$\phi 21.1 \pm 0.3$</td> <td>$\phi 32.2 \pm 0.2$</td> <td>15.5 ± 0.2</td> <td>83.5 ± 0.5</td> </tr> <tr> <td>PE</td> <td>$\phi 38.6 \pm 0.3$</td> <td>81.65 ± 0.3</td> <td>$\phi 20.9 \pm 0.3$</td> <td>$\phi 32.2 \pm 0.2$</td> <td>15.5 ± 0.2</td> <td>83.65 ± 0.5</td> </tr> </tbody> </table>		location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	total H (mm)	60ml	PP	$\phi 39.0 \pm 0.3$	81.5 ± 0.3	$\phi 21.1 \pm 0.3$	$\phi 32.2 \pm 0.2$	15.5 ± 0.2	83.5 ± 0.5	PE	$\phi 38.6 \pm 0.3$	81.65 ± 0.3	$\phi 20.9 \pm 0.3$	$\phi 32.2 \pm 0.2$	15.5 ± 0.2
location name	material	O. D. D1 (mm)	Height H1 (mm)	I. D. D2 (mm)	Cap O.D D3 (mm)	Cap H H2 (mm)	total H (mm)																	
60ml	PP	$\phi 39.0 \pm 0.3$	81.5 ± 0.3	$\phi 21.1 \pm 0.3$	$\phi 32.2 \pm 0.2$	15.5 ± 0.2	83.5 ± 0.5																	
	PE	$\phi 38.6 \pm 0.3$	81.65 ± 0.3	$\phi 20.9 \pm 0.3$	$\phi 32.2 \pm 0.2$	15.5 ± 0.2	83.65 ± 0.5																	
Colour	Amber																							
Nominal Capacity	60ml																							
Material Characteristics	bottle:High Density Polyethylene(HDPE)																							
	cap:polypropylene (PP)																							
Performance parameters	Negative pressure:-70kpa~-95kpa																							
	Ethanol loss rate: $\leq 0.5\%$																							
	water loss rate: $\leq 0.2\%$																							
Purity Grade	free of RNase and DNase																							
	free of pyrogen																							
	free of PCR inhibitors																							
	free of human DNA																							
Shelf - Life	60 months at controlled storage conditions																							
Packaging	pcs per case:800																							
Storage	dry condition:40%-60% rel. humidity																							
Approved by PakGent quality department																								