


Technical Datasheet

Article - No:	TDS-NMPB250PA	Revision: A0																						
Remarks	single use only																							
Specification																								
Product name	NMPB250PA,250ml Packaging Bottle,PP,Amber																							
Production conditions	Class 100,000 clean room production																							
	No release agent is used																							
Classification	CE/IVDD																							
Dimensions																								
	<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th>location name</th> <th>material</th> <th>O.D. D1 (mm)</th> <th>Height H1 (mm)</th> <th>I.D. D2 (mm)</th> <th>Cap O.D. D3 (mm)</th> <th>Cap H H2 (mm)</th> <th>Total H (mm)</th> </tr> </thead> <tbody> <tr> <td rowspan="2">250ml</td> <td>PP</td> <td>ø60.8</td> <td>130.0</td> <td>ø17.5</td> <td>ø28.4</td> <td>15.1</td> <td>131.2</td> </tr> <tr> <td>PE</td> <td>ø60.8</td> <td>130.0</td> <td>ø17.5</td> <td>ø28.4</td> <td>15.1</td> <td>131.2</td> </tr> </tbody> </table>		location name	material	O.D. D1 (mm)	Height H1 (mm)	I.D. D2 (mm)	Cap O.D. D3 (mm)	Cap H H2 (mm)	Total H (mm)	250ml	PP	ø60.8	130.0	ø17.5	ø28.4	15.1	131.2	PE	ø60.8	130.0	ø17.5	ø28.4	15.1
location name	material	O.D. D1 (mm)	Height H1 (mm)	I.D. D2 (mm)	Cap O.D. D3 (mm)	Cap H H2 (mm)	Total H (mm)																	
250ml	PP	ø60.8	130.0	ø17.5	ø28.4	15.1	131.2																	
	PE	ø60.8	130.0	ø17.5	ø28.4	15.1	131.2																	
Colour	Amber																							
Nominal Capacity	250ml																							
Material Characteristics	bottle:polypropylene (PP)																							
	cap:High Density Polyethylene(HDPE)																							
Performance parameters	Negative pressure:-70kpa~-95kpa																							
	Ethanol loss rate: ≤0.5%																							
	water loss rate: ≤0.2%																							
Purity Grade	free of RNase and DNase																							
	free of pyrogen																							
	free of PCR inhibitors																							
	free of human DNA																							
	free of endotoxins																							
Shelf - Life	60 months at controlled storage conditions																							
Packaging	pcs per case:250																							
Storage	dry condition:40%-60% rel. humidity																							
Approved by PakGent quality department																								