


## Technical Datasheet

<b>Article - No:</b>	TDS-CL-FT075V	Revision: A0
<b>Remarks</b>	single use only	
Specification		
<b>Product name</b>	CL-FT075V,75cm <sup>2</sup> ,Cell Culture Flask,TC Treated,Vent Cap	
<b>Production conditions</b>	Class 10,000 clean room production	
	No release agent is used	
<b>Classification</b>	CE/IVD	
<b>Dimensions</b>		
<b>Colour</b>	Clear	
<b>Graduation</b>	moulded rings	
<b>Nominal Capacity</b>	75cm <sup>2</sup>	
<b>Material Characteristics</b>	bottle:polystyrene (PS)	
	cap:High Density Polyethylene(HDPE)	
<b>Sterile</b>	Yes	
<b>Product Characteristics</b>	biologically inert	
	no heavy metals according to EC/94/62	
<b>Purity Grade</b>	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
<b>Shelf - Life</b>	60 months at controlled storage conditions	
<b>Packaging</b>	pcs per bag: 5	
	bag per case: 20	
<b>Storage</b>	Storage Temperature:20°C-25°C,humidity40%-60%	
	Keep clean,ventilated and out of direct sunlight	
<b>Approved by PakGent quality department</b>		