


Technical Datasheet

Article - No:	TDS-CL-F225V	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-F225V, 225cm ² , Cell Culture Flask, Vent cap	
Production conditions	Class 10,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	225cm ²	
Material Characteristics	bottle: polystyrene (PS)	
	cap: High Density Polyethylene (HDPE)	
Sterile	Yes	
Product Characteristics	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	60 months at controlled storage conditions	
Packaging	pcs per bag: 5	
	bag per case: 5	
Storage	Storage Temperature: 20°C-25°C, humidity 40%-60%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		