


Technical Datasheet

Article - No:	TDS-PCRP-20NS-W	Revision: A0																						
Remarks	single use only																							
Specification																								
Product name	PCRP-20NS-W,0.2ml PCR 96-well plate																							
Production conditions	Class 100,000 clean room production																							
	No release agent is used																							
Classification	CE/IVD																							
Dimensions	 <table border="1" style="margin-left: auto; margin-right: auto; border-collapse: collapse; text-align: center;"> <thead> <tr> <th>A</th><th>B</th><th>C</th><th>D</th><th>E</th><th>F</th><th>G</th><th>H</th><th>I</th><th>J</th><th>K</th> </tr> </thead> <tbody> <tr> <td>99.19</td><td>84.94</td><td>125.27</td><td>63.12</td><td>19.47</td><td>∅7.12</td><td>21.12</td><td>17.50°</td><td>1.02</td><td>R1.97</td><td>C7.9</td> </tr> </tbody> </table>		A	B	C	D	E	F	G	H	I	J	K	99.19	84.94	125.27	63.12	19.47	∅7.12	21.12	17.50°	1.02	R1.97	C7.9
A	B	C	D	E	F	G	H	I	J	K														
99.19	84.94	125.27	63.12	19.47	∅7.12	21.12	17.50°	1.02	R1.97	C7.9														
Colour	White																							
Nominal Capacity	0.2ml																							
With skirt	non-skirt																							
Material Characteristics	polypropylene (PP), premium surface																							
Product Characteristics	biologically inert																							
	no heavy metals according to EC/94/62																							
Purity Grade	free of RNase and DNase																							
	free of pyrogen																							
	free of PCR inhibitors																							
	free of human DNA																							
Shelf - Life	36 months at controlled storage conditions																							
	36 months at controlled storage conditions																							
Packaging	pcs per box: 10																							
	box per case: 10																							
Storage	dry condition:40%-60% rel. humidity																							
Approved by PakGent quality department																								